SYMBOLS DESCRIPTION

Full explanation of symbols used on Kerr packaging located at:
<http://www.kerrdental.com/symbols-glossary>

REF	Catalogue number
Manufacture	
LOT	Batch Code
CE 2797	European conformity
Consult instructions for use	
Caution	
Authorized representative in the European Community	
Use by date	
Do not re-use	
Keep away from sunlight	
MD	1. Medical Device
28°C	2. Temperature Storage Limitations: Low 2°C/36°F and High 28°C/82°F

CONTENTS

en - ENGLISH

OptiShade™

Nanohybrid Universal Restorative Composite
ISO 4049: Polymer-Based Restorative Material, Class 2/Group 1

GENERAL INFORMATION

OptiShade is a light-cured, esthetic, resin-based nanohybrid dental restorative designed for the placement of direct restorations. OptiShade is designed to be used after placement of a dental adhesive, such as a product from the Kerr OptiBond® family. OptiShade is available in Unidose composite compules and syringe delivery systems.

This product is intended to be used by qualified dental professionals on a general patient population.

PRODUCT COMPOSITION

100% aluminum has a radio-opacity equivalent to that of dentin and 200% aluminum is equivalent to enamel.

The restorative is available in three composite shades (Light, Medium, Dark) and two special shades (Bleach White and Universal Opaque). All shades are radiopaque with about 250% aluminum.

The fillers comprise selection of chemically infused mixed oxides, pre-polymerized filler, barium glass filler, silica, and yttrium trifluoride (about 61.5% weight or 65.1% by volume for the Light, Medium, Dark, and Universal Opaque shades; about 78.5% weight or 60% by volume for the Bleach White shade). The inorganic fillers are up to 75% weight or 55% by volume. The Light, Medium, Dark, and Universal Opaque shades have a range particle size of 40 nm to 6 microns. The Bleach White shade has a range particle size of 40 nm to 30 microns. OptiShade contains Bis-EMA, Bis-GMA, and TEGDMA resins.

INDICATIONS FOR USE
OptiShade, a nanohybrid dental restorative, is indicated for direct placement in all cavity classes in anterior and posterior teeth. Additional indications include: repair of porcelain/composite restorations, and core build-ups.

CONTRAINDICATIONS

None known.

WARNINGS

1. Do not use this product on patients with known methacrylate allergies. Uncured methacrylic resin may cause contact dermatitis and damage to the pulp.
2. Avoid contact with skin, eyes, and soft tissue. Wash thoroughly with water after contact.
3. Get medical attention for eyes.
4. Do not take internally.

PRECAUTIONS

Kerr recommends the following precautions when using OptiShade or any other composite material:

- Provide the patient with protective eyewear.
- Use a rubber dam to isolate the oral cavity.
- Avoid fluid contamination during placement.
- Use of protective gloves and a no-touch clinical technique is strongly recommended. If the composite is touched with commonly used gloves, acrylate may penetrate the glove. Dispose of glove immediately and wash hands.
- Protect the material from ambient light that may initiate polymerization.
- The Unidose tips are designed for SINGLE PATIENT USE ONLY to prevent cross-contamination between patients.

Kerr, OptiShade, and Unidose are trademarks owned by Kerr Corporation.

ADVERSE EVENTS

If a serious incident occurs with this medical device, report it to the manufacturer and to the competent authority for that country in which the user and/or patient is established.

STEEP-BY-STEP INSTRUCTIONS

1. Wear personal protective equipment when using and handling this product.
2. Do not use this product near eugenol-containing products.
3. Clean the tooth with non-fluoridated pumice and water to remove any debris or stain that may interfere with the shade selection process.

Select the desired shade based on the chart listed below. A composite mock-up or VITA®-matched shade guide may be used to check the color of the tooth to be restored and any adjacent teeth. The shade of the restorative material is affected by the thickness of the composite shade selection and mock-ups should match the thickness of the tooth to be restored.

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